

Anti-SKP2 (4D8) Mouse antibody



Catalog Number: 168392

Product name

Anti-SKP2 (4D8) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human

Applications

WB

Antibody description

SKP2 (4D8) Mouse monoclonal antibody

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000


Preparation

Antigen: Purified recombinant fragment of SKP2 (aa1-130) expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Skp2 in HeLa cell lysates using Skp2 mouse mAb (1:2000 diluted). Predicted band size: 48KDa. Observed band size: 48KDa.

Clonality

Monoclonal