

Anti-Hsp70 (C-term) (2A11) Mouse antibody



Catalog Number: 168401

Product name

WB, IP

Anti-Hsp70 (C-term) (2A11) Mouse antibody

Dilutions

Specificity

WB: 1/1000

Human

Antibody description

Hsp70 (C-term) (2A11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Hsp70(C-terminus) protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine(pH 7.4 150 mM NaCl)with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG2b

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysates using Hsp70 (C-terminus) mouse mAb.



Western blot detection of Hsp70(C-terminus) in HeLa and HsHeLa cell lysates using Hsp70(C-terminus) mouse mAb (1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.