

Anti-PRMT1 (4E9) Mouse antibody



Catalog Number: 168402

Product name

Anti-PRMT1 (4E9) Mouse antibody

Specificity

Human

Antibody description

PRMT1 (4E9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human PRMT1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of PRMT1 in HeLa, A549 and SW480 cell lysates using PRMT1 mouse mAb (1:1000 diluted). Predicted band size: 42 kDa. Observed band size: 42 kDa.