

Anti-AIF (8H1) Mouse antibody



Catalog Number: 168405

Product name

Anti-AIF (8H1) Mouse antibody

Specificity

Human

Antibody description

AIF (8H1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human AIF protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/200

Validations



Western blot analysis of extracts from HeLa, Ramos, HepG2, MCF7 and Jurkat cell lysates using AIF mouse mAb (1:1000 diluted). Predicted band size: 67KDa. Observed band size: 67KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-AIF mouse mAb (dilution 1:200).