

# Anti-CTNNBL1 Mouse antibody



Catalog Number: 168412

## Product name

Anti-CTNNBL1 Mouse antibody

WB: 1/500 - 1/2000

FC: 1/200 - 1/400

## Specificity

Human

ELISA: 1/10000

## Antibody description

CTNNBL1 Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human CTNNBL1 (AA: 390-557) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


Mouse IgG2a

## Applications

WB, FC


## Dilutions

## Validations



Flow cytometric analysis of HeLa cells using CTNNBL1 mouse mAb (green) and negative

control (red).



Western blot detection of CTNNBL1 in A549, MCF7 and HeLa cell lysates using CTNNBL1 mouse mAb (1:500 diluted). Predicted band size: 65 kDa. Observed band size: 65 kDa.