

Anti-Glutathione S Transferase mu 1 (1H4) Mouse antibody



Catalog Number: 168416

Product name

Mouse IgG1

Anti-Glutathione S Transferase mu 1 (1H4) Mouse antibody

Applications

WB, IHC-P, FC

Specificity

Human, Monkey

Dilutions

WB: 1/500 - 1/2000

Antibody description

Glutathione S Transferase mu 1 (1H4) Mouse monoclonal antibody

IHC: 1/200 - 1/1000

FC: 1/200 - 1/400

Preparation

Antigen: Purified recombinant fragment of human GSTM1 expressed in E. Coli.

ELISA: 1/10000

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Glutathione S Transferase mu in MCF7, Jurkat and HeLa cell lysates using Glutathione S Transferase mu mouse mAb (dilution 1:1000). Predicted band size: 26 Kda. Observed band size: 26KDa.