

Anti-IRAK4 (2H9) Mouse antibody



Catalog Number: 168423

Product name

Anti-IRAK4 (2H9) Mouse antibody

Specificity

Human, Mouse, Monkey

Antibody description

IRAK4 (2H9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human IRAK4 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, FC

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

Validations



Western blot detection of IRAK4 in MCF7 cell lysates using IRAK4 mouse mAb (dilution 1:500). Predicted band size: 52 Kda. Observed band size: 55KDa.