

Anti-RPA32/RPA2 (3E7) Mouse antibody



Catalog Number: 168447

Product name

Anti-RPA32/RPA2 (3E7) Mouse antibody

Specificity

Human

Antibody description

RPA32/RPA2 (3E7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human RPA32/RPA2 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, IP, ICC/IF

Dilutions


WB: 1/2000

;ICC: 1/200


Validations



Immunocytochemistry staining of HeLa cells fixed with -20u2103 Methanol and using anti-RPA32/RPA2 antibody (dilution 1:200).



Immunoprecipitation analysis of HeLa cell lysates using RPA32/RPA2 mouse mAb.



Western blot detection of RPA32/RPA2 in MCF7, HeLa and HUVEC cell lysates using RPA32/RPA2 mouse mAb (1:2000 diluted). Predicted band size: 32 kDa. Observed band size: 32 kDa. Exposure time: 20s.