

Anti-Cytochrome P450 17A1/CYP17A1 (6E1) Mouse antibody



Catalog Number: 168456

Product name

Anti-Cytochrome P450 17A1/CYP17A1 (6E1)
Mouse antibody

Specificity

Human

Antibody description

Cytochrome P450 17A1/CYP17A1 (6E1) Mouse
monoclonal antibody

Preparation

Antigen: Purified recombinant human Cytochrome
P450 17A1 protein fragments expressed in E.coli.

Formulation

Ascites containing with 0.03% Proclin300 and 50%
glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Cytochrome P450 17A1
in HeLa, HepG2 and A549 cell lysates using
Cytochrome P450 17A1 mouse mAb (1:1000
diluted). Predicted band size: 60KDa. Observed
band size: 60KDa.