

Anti-MSH6 (3H9) Mouse antibody



Catalog Number: 168460

Product name

Anti-MSH6 (3H9) Mouse antibody

Dilutions

WB: 1/500 - 1/2000

Specificity

ELISA: 1/10000

Human

Antibody description

MSH6 (3H9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human MSH6 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB

Validations



Western blot analysis using MSH6 mouse mAb against MCF-7 (1), HEK293 (2), and HCT116 (3)



cell lysate.

Western blot detection of MSH6 in K562, 3T3 and HeLa cell lysates using MSH6 mouse mAb (1:500 diluted). Predicted band size: 160 kDa. Observed band size: 160 kDa.