

# Anti-JNK1 (1H8) Mouse antibody



Catalog Number: 168465

## Product name

Anti-JNK1 (1H8) Mouse antibody

## Specificity

Human, Mouse

## Antibody description

JNK1 (1H8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human MAPK8 expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, FC


## Dilutions

WB: 1/500 - 1/2000


FC: 1/200 - 1/400

ELISA: 1/10000

## Validations



Flow cytometric analysis of HeLa cells using MAPK8 mouse mAb (green) and negative control



Western blot analysis of extracts from Raji, HCT116, Jurkat, MCF7 and A549 cell lysates using JNK1 mouse mAb (1:1000 diluted). Predicted band size: 48 kDa. Observed band size: 46, 54 kDa.