

Anti-EEF2/Elongation factor 2 (4B3) Mouse antibody



Catalog Number: 168473

Product name

Anti-EEF2/Elongation factor 2 (4B3) Mouse antibody

Specificity

Human

Antibody description

EEF2/Elongation factor 2 (4B3) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human eEF2 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, ICC/IF, ChIP

Dilutions


WB: 1/5000

;ICC: 1/200

Validations



Western blot detection of eEF2 in HL-60, Jurkat, SH-SY-5Y, U2OS and HeLa cell lysates using eEF2 mouse mAb (1:5000 diluted). Predicted band size: 95 kDa. Observed band size: 95 kDa. Exposure time: 15s.



Immunocytochemistry staining of HeLa cells fixed with -20u2103 Methanol and using anti-eEF2 mouse mAb (dilution 1:200).