

Anti-ALDH2 Mouse antibody



Catalog Number: 168478

Product name

Anti-ALDH2 Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

ALDH2 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human ALDH2 (AA: 317-517) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

Validations



Western blot detection of ALDH2 in A549, MCF7, C6 and HeLa cell lysates using ALDH2 mouse mAb (1:3000 diluted). Predicted band size: 56 kDa. Observed band size: 56 kDa.