

Anti-KAT5/Tip60 Mouse antibody



Catalog Number: 168484

Product name

Mouse IgG1

Anti-KAT5/Tip60 Mouse antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human

Dilutions

Antibody description

WB: 1/500 - 1/2000

KAT5/Tip60 Mouse monoclonal antibody

FC: 1/200 - 1/400

Preparation

ELISA: 1/10000

Antigen: Purified recombinant fragment of human TIP60 (AA: 18-208) expressed in E. Coli.

Validations

Formulation

Purified antibody in PBS with 0.05% sodium azide


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Tip60 in MCF7, C6 and HeLa cell lysates using Tip60 mouse mAb (1:2000 diluted). Predicted band size: 60KDa. Observed band size: 60KDa.