

Anti-LYRIC/AEG1 (2E4) Mouse antibody



Catalog Number: 168492

Product name

Anti-LYRIC/AEG1 (2E4) Mouse antibody

Specificity

Human

Antibody description

LYRIC/AEG1 (2E4) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human Metadherin expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, FC

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000


FC: 1/200 - 1/400

ELISA: 1/10000

Validations




Immunohistochemical analysis of paraffin-embedded human Liver tissues using Metadherin



mouse mAb

Western blot detection of Lyric/Metadherin in Jurkat, A549, MCF7 and HeLa cell lysates using Lyric/Metadherin mouse mAb (1:5000 diluted). Predicted band size: 64KDa. Observed band size: 75KDa.



Flow cytometric analysis of HeLa cells using Metadherin mouse mAb (green) and negative control (purple).