

# Anti-SQSTM1/p62 (3F8) Mouse antibody



Catalog Number: 168499

## Product name

Anti-SQSTM1/p62 (3F8) Mouse antibody

## Specificity

Human

## Antibody description

SQSTM1/p62 (3F8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human SQSTM1 (AA: 232-356) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, FC

## Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

## Validations



Western blot analysis of extracts from Jurkat, A549, MCF7, C6 and HeLa cell lysates using SQSTM1/p62 mouse mAb (1:1000 diluted). Predicted band size: 62 kDa. Observed band size: 62 kDa.