

# Anti-HSP40 Mouse antibody



Catalog Number: 168506

## Product name

Anti-HSP40 Mouse antibody

## Specificity

Human

## Antibody description

HSP40 Mouse monoclonal antibody

## Preparation

Antigen: Recombinant Protein of Hsp40

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at -20°C

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

1:1000

## Validations



Western blot detection of HSP40 in Jurkat cell lysates using HSP40 Mouse mAb(1:1000 diluted). Predicted band size:40KDa. Observed band size:40KDa.