

Anti-KAT13D/CLOCK Mouse antibody



Catalog Number: 168507

Product name

Anti-KAT13D/CLOCK Mouse antibody

Specificity

Human

Antibody description

KAT13D/CLOCK Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CLOCK expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ICC/IF

Dilutions


ICC: 1/200 - 1/1000

ELISA: 1/10000

Validations



Immunofluorescence analysis of HeLa cells using CLOCK mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Western blot analysis using CLOCK mouse mAb against CLOCK(AA: 200-465)-hIgGfc transfected HEK293 cell lysate.