

# Anti-STAT3(1E3) Mouse antibody

Catalog Number: 168511

## Product name

Anti-STAT3(1E3) Mouse antibody

## Specificity

Human, Rat, Mouse

## Antibody description

STAT3(1E3) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant Protein of STAT3

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

-20°C/1 year

## Clonality

Monoclonal

## Ig Type

Unknown

## Applications


WB, IHC-P

## Dilutions


1:1000-2000

1:100-200


## Validations



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma using STAT3 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using STAT3 Mouse mAb diluted at 1:200.



Western blot analysis of 1) HeLa, 2) 3T3, 3) PC12 with STAT3 Mouse mAb diluted at 1:2,000.