

Anti-GFAP Mouse antibody



Catalog Number: 168520

Product name

Anti-GFAP Mouse antibody

Specificity

Human

Antibody description

GFAP Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human GFAP expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000


ICC: 1/200 - 1/1000

ELISA: 1/10000

Validations




Immunofluorescence analysis of paraffin-embedded lobe of brain tissues using GFAP mouse



mAb (green).

Western blot analysis using GFAP mouse mAb against A431 (1), SK-N-SH (2) and PC12 (3) cell



lysate.

Immunohistochemical analysis of paraffin-embedded brain tissues using GFAP mouse mAb with DAB staining