

Anti-Dynein Light Chain (4D12) Mouse antibody



Catalog Number: 168527

Product name

Mouse IgG1

Anti-Dynein Light Chain (4D12) Mouse antibody

Applications

Specificity

WB, FC

Human

Dilutions

Antibody description

WB: 1/500 - 1/2000

Dynein Light Chain (4D12) Mouse monoclonal antibody

FC: 1/200 - 1/400

ELISA: 1/10000

Preparation

Antigen: Purified recombinant fragment of human DNAL4 (AA: 1-105) expressed in E. Coli.

Validations

Formulation

Purified antibody in PBS with 0.05% sodium azide


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot analysis of extracts from C6 cell lysates using Dynein light chain mouse mAb (1:500 diluted). Predicted band size: 12 kDa. Observed band size: 12 kDa.