

# Anti-Bid (3C5) Mouse antibody



Catalog Number: 168529

## Product name

Anti-Bid (3C5) Mouse antibody

IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

## Specificity

Human

FC: 1/200 - 1/400

ELISA: 1/10000

## Antibody description

Bid (3C5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human BID expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, IHC-P, ICC/IF, FC


## Dilutions

WB: 1/500 - 1/2000

## Validations



Immunocytochemistry staining of HeLa cells using Bid mouse mAb (dilution 1:100).



Western blot detection of Bid in 293T, Jurkat and A431 cell lysates using Bid mouse mAb (1:300 diluted). Predicted band size: 22KDa. Observed band size: 22KDa.