

Anti-BRAF (3E8) Mouse antibody



Catalog Number: 168533

Product name

Anti-BRAF (3E8) Mouse antibody

Specificity

Human, Mouse

Antibody description

BRAF (3E8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human BRAF expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000


ICC: 1/200 - 1/1000

ELISA: 1/10000


Validations



Immunohistochemical analysis of paraffin-embedded human testis tissues using BRAF mouse



Western blot detection of B Raf in A549, Ramos, Jurkat, THP-1 and Raji cell lysates using B Raf mouse mAb (1:1000 diluted). Predicted band size: 87KDa. Observed band size: 87KDa.



Confocal immunofluorescence analysis of MCF-7 (left) and HepG2 (right) cells using BRAF mouse mAb (green). Red