

Anti-PYK2 (4B4) Mouse antibody



Catalog Number: 168534

Product name

Anti-PYK2 (4B4) Mouse antibody

WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

Specificity


Human

ELISA: 1/10000

Antibody description

PYK2 (4B4) Mouse monoclonal antibody

Validations



Preparation

Antigen: Purified recombinant fragment of PYK2 (aa815-997) expressed in E. Coli.

Western blot detection of Pyk2 in A549, Ramos, Jurkat, THP-1 and Raji cell lysates using Pyk2 mouse mAb (1:1000 diluted). Predicted band size: 116KDa. Observed band size: 116KDa.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2a

Applications

WB, IHC-P

Dilutions



Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using PYK2 mouse mAb.