

# Anti-Moesin (2C12) Mouse antibody



Catalog Number: 168540

## Product name

Anti-Moesin (2C12) Mouse antibody

## Specificity

Human, Monkey

## Antibody description

Moesin (2C12) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human MSN expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, FC

## Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

## Validations



Western blot detection of Moesin in COS7, Jurkat, HaCat and HeLa cell lysates using Moesin mouse mAb (dilution 1:1000). Predicted band size: 67.8KDa. Observed band size: 67.8KDa.