

Anti-GRP78 BiP (3B4) Mouse antibody



Catalog Number: 168541

Product name

WB: 1/500 - 1/2000

Anti-GRP78 BiP (3B4) Mouse antibody

IHC: 1/200 - 1/1000

Specificity

ELISA: 1/10000

Human

Antibody description

GRP78 BiP (3B4) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human HSPA5 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications


WB, IHC-P

Dilutions

Validations



Western blot analysis of extracts from Jurkat, A549, MCF7 and C6 cell lysates using BiP/GRP78 mouse mAb (1:1000 diluted). Predicted band size: 72KDa. Observed band size: 72KDa.



Immunohistochemical analysis of paraffin-embedded human Thyroid tissues using HSPA5 mouse mAb