

Anti-KBTBD8 Mouse antibody



Catalog Number: 168550

Product name

Anti-KBTBD8 Mouse antibody

Specificity

Human

Antibody description

KBTBD8 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human KBTBD8 (AA: 264-464) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF, FC


Dilutions

WB: 1/500 - 1/2000

FC: 1/200 - 1/400

ELISA: 1/10000

Validations



Western blot detection of KBTBD8 in Jurkat, A549, MCF7, C6 and HeLa cell lysates using KBTBD8 mouse mAb (1:1000 diluted). Predicted band size: 69KDa. Observed band size: 69KDa.