

# Anti-PKR Rabbit antibody



Catalog Number: 168565

## Product name

Anti-PKR Rabbit antibody

## Specificity

Human

## Antibody description

PKR Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human PKR

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

1:2000

1:20

1:20


1:20

1:20

## Validations



Immunocytochemistry of PKR(green) in HeLa cells using PKR Rabbit mAb at dilution 1/200, and



Western blot detection of PKR in K562, HeLa cell lysates using PKR Rabbit mAb(1:1000 diluted). Predicted band size: 62kDa. Observed band size: 74kDa.