

# Anti-GAP43 Rabbit antibody



Catalog Number: 168568

## Product name

Anti-GAP43 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GAP43 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human GAP43

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

1:5000-1:20000


1:50

1:50

1:20

1:20

## Validations



Western blot detection of GAP43 in Rat Brain lysates using GAP43 Rabbit mAb(1:1000 diluted). Predicted band size:25kDa. Observed band size:46kDa.