

Anti-Histone H4 Rabbit antibody



Catalog Number: 168570

Product name

Rabbit IgG

Anti-Histone H4 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, ChIP

Human, Mouse, Rat

Dilutions

Antibody description

1:1000

Histone H4 Rabbit monoclonal antibody

1:200

Preparation

1:20


Antigen: A synthetic peptide of human Histone H4

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Histone H4 in HeLa cell lysates using Histone H4 Rabbit mAb(1:1000 diluted).Predicted band size:11kDa.Observed band size:11kDa.

Clonality

Monoclonal

Ig Type