

Anti-FTO Rabbit antibody



Catalog Number: 168575

Product name

Anti-FTO Rabbit antibody

Specificity

Human

Antibody description

FTO Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human FTO

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

1:1000

1:20

Validations



Western blot detection of FTO in K562, HeLa cell lysates using FTO Rabbit mAb(1:1000 diluted). Predicted band size: 58kDa. Observed band size: 58kDa.