

# Anti-Adenosine A1 Receptor Rabbit antibody



Catalog Number: 168578

## Product name

Anti-Adenosine A1 Receptor Rabbit antibody

## Specificity

Human

## Antibody description

Adenosine A1 Receptor Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Adenosine A1 Receptor

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

1:1000


1:50

1:20

1:20-1:100

1:20

## Validations



Western blot detection of Adenosine A1 Receptor in K562,C6,3T3,HeLa cell lysates using Adenosine A1 Receptor Rabbit mAb(1:1000 diluted).Predicted band size:37kDa.Observed band size:37kDa.