

# Anti-PP-X Rabbit antibody



Catalog Number: 168588

## Product name

Anti-PP-X Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PP-X Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human PPP4C

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

1:1000

1:20

1:50

## Validations



Western blot detection of PP-X in K562, C6, 3T3, HeLa cell lysates using PP-X Rabbit mAb(1:1000 diluted). Predicted band size: 35kDa. Observed band size: 35kDa.