

# Anti-LAT Rabbit antibody



Catalog Number: 168608

## Product name

Anti-LAT Rabbit antibody

## Specificity

Human

## Antibody description

LAT Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human LAT

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000

IHC: 1/20

FC: 1/20

IP: 1/20

## Validations



Western blot detection of LAT in U2OS cell lysates using LAT Rabbit mAb(1:500 diluted).Predicted band size:28kDa.Observed band size:36kDa.