

# Anti-p18 INK4c/CDKN2C Rabbit antibody



Catalog Number: 168613

## Product name

Rabbit IgG

Anti-p18 INK4c/CDKN2C Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000-1/2000

p18 INK4c/CDKN2C Rabbit monoclonal antibody

IHC: 1/20

## Preparation

ICC/IF: 1/20


Antigen: Recombinant protein of human p18 INK4c

IP: 1/20

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of p18 INK4c/CDKN2C in C6 cell lysates using p18 INK4c/CDKN2C Rabbit mAb(1:500 diluted).Predicted band size:18kDa.Observed band size:18kDa.

## Clonality

Monoclonal

## Ig Type