

# Anti-MAS1L Rabbit antibody



Catalog Number: 168615

## Product name

Anti-MAS1L Rabbit antibody

## Specificity

Human

## Antibody description

MAS1L Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human MAS1L

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

1:2000-1:10000

1:20

1:20-1:100

## Validations



Western blot detection of MAS1L in K562,Rat Brain,C6,3T3,Hela cell lysates using MAS1L Rabbit mAb(1:1000 diluted).Predicted band size:42kDa.Observed band size:42kDa.