

# Anti-CD55 Rabbit antibody



Catalog Number: 168623

## Product name

Anti-CD55 Rabbit antibody

## Specificity

Human

## Antibody description

CD55 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human CD55

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/2000-1/10000

## Validations



Western blot detection of CD55 in HeLa cell lysates using CD55 Rabbit mAb(1:1000 diluted).Predicted band size:41kDa.Observed band size:78kDa.