

Anti-AMPK alpha 2 Rabbit antibody



Catalog Number: 168625

Product name

Anti-AMPK alpha 2 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Rat

Applications

WB, IP

Antibody description

AMPK alpha 2 Rabbit monoclonal antibody

Dilutions

1:1000

Preparation

Antigen: Recombinant protein of human AMPK alpha 2

1:20

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of AMPK alpha 2 in C6, HeLa cell lysates using AMPK alpha 2 Rabbit mAb(1:1000 diluted). Predicted band size: 62kDa. Observed band size: 62kDa.