

# Anti-UFD1L Rabbit antibody



Catalog Number: 168626

## Product name

Anti-UFD1L Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

UFD1L Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human UFD1L

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000-1/5000

## Validations



Western blot detection of UFD1L in HeLa cell lysates using UFD1L Rabbit mAb(1:1000 diluted).Predicted band size:35kDa.Observed band size:40kDa.