

# Anti-TdT Rabbit antibody



Catalog Number: 168637

## Product name

Anti-TdT Rabbit antibody

## Specificity

Human

## Antibody description

TdT Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human TdT

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/2000

## Validations



Western blot detection of TdT in Jurkat cell lysates using TdT Rabbit mAb(1:1000 diluted).Predicted band size:59kDa.Observed band size:59kDa.