

# Anti-UHRF1 Rabbit antibody



Catalog Number: 168638

## Product name

Anti-UHRF1 Rabbit antibody

## Specificity

Human, Rat

## Antibody description

UHRF1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human UHRF1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/50

FC: 1/20

IP: 1/20

## Validations



Western blot detection of UHRF1 in K562, HeLa cell lysates using UHRF1 Rabbit mAb(1:1000 diluted). Predicted band size: 89kDa. Observed band size: 97kDa.