

# Anti-CD31 Rabbit antibody



Catalog Number: 168639

## Product name

Anti-CD31 Rabbit antibody

## Specificity

Human

## Antibody description

CD31 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human CD31

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/20

## Validations



Western blot detection of CD31 in Jurkat,3T3 cell lysates using CD31 Rabbit mAb(1:1000 diluted).Predicted band size:83kDa.Observed band size:130kDa.