

Anti-SOX9 Rabbit antibody



Catalog Number: 168640

Product name

Anti-SOX9 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SOX9 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human SOX9

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/2000

ICC/IF: 1/20

FC: 1/20

Validations



Western blot detection of SOX9 in C6,3T3,HeLa cell lysates using SOX9 Rabbit mAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:70kDa.