

# Anti-HNF-4-alpha Rabbit antibody



Catalog Number: 168655

## Product name

Rabbit IgG

Anti-HNF-4-alpha Rabbit antibody

## Applications

## Specificity

WB, IP, ChIP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

HNF-4-alpha Rabbit monoclonal antibody

IP: 1/20

## Preparation

ChIP: 1/20

Antigen: Recombinant protein of human HNF-4-alpha

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of HNF-4-alpha in K562,C6,HeLa cell lysates using HNF-4-alpha Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.