

Anti-INPP4B Rabbit antibody



Catalog Number: 168657

Product name

Anti-INPP4B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

INPP4B Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human INPP4B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/20

Validations



Western blot detection of INPP4B in Rat Brain, C6 cell lysates using INPP4B Rabbit mAb(1:1000 diluted). Predicted band size: 105kDa. Observed band size: 105kDa.