

Anti-TNFAIP8 Rabbit antibody



Catalog Number: 168661

Product name

Anti-TNFAIP8 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP

Antibody description

TNFAIP8 Rabbit monoclonal antibody

Dilutions

1:1000

1:20

Preparation

Antigen: A synthetic peptide of human TNFAIP8

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of TNFAIP8 in CEM, K562 cell lysates using TNFAIP8 Rabbit mAb(1:1000 diluted). Predicted band size: 23kDa. Observed band size: 21-23kDa.