

Anti-Thrombomodulin Rabbit antibody



Catalog Number: 168678

Product name

Anti-Thrombomodulin Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse, Rat

Applications

WB, FC

Antibody description

Thrombomodulin Rabbit monoclonal antibody

Dilutions

1:1000

Preparation


Antigen: Recombinant protein of mouse
Thrombomodulin

1:20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Thrombomodulin in 3T3 cell lysates using Thrombomodulin Rabbit mAb(1:1000 diluted).Predicted band size:62kDa.Observed band size:95kDa.

Clonality

Monoclonal