

Anti-DNAJA4 Rabbit antibody



Catalog Number: 168689

Product name

Anti-DNAJA4 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

DNAJA4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human DNAJA4

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/1000-1/2000

FC: 1/50

IP: 1/20

Validations



Western blot detection of DNAJA4 in C6,Rat Brain lysates using DNAJA4 Rabbit mAb(1:500 diluted).Predicted band size:45kDa.Observed band size:45kDa.