

Anti-JNK1/2/3 Rabbit antibody



Catalog Number: 168702

Product name

Rabbit IgG

Anti-JNK1/2/3 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

1:1000

JNK1/2/3 Rabbit monoclonal antibody

1:200

Preparation

1:20


Antigen: Recombinant protein of human JNK1

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of JNK1/2/3 in Jurkat,Rat Brain,C6,CHO-K1,Hela cell lysates using JNK1/2/3 Rabbit mAb(1:1000 diluted).Predicted band size:48,53 kDa.Observed band size:46,54kDa.

Clonality

Monoclonal

Ig Type