

Anti-PCPTP1/PTPRQ Rabbit antibody



Catalog Number: 168706

Product name

Anti-PCPTP1/PTPRQ Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

PCPTP1/PTPRQ Rabbit monoclonal antibody

Dilutions

1:1000

Preparation

Antigen: A synthetic peptide of human PCPTP1

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of PCPTP1/PTPRQ in Rat Brain lysates using PCPTP1/PTPRQ Rabbit mAb(1:1000 diluted).Predicted band size:74kDa.Observed band size:74kDa.